





Guide to choosing the Inflatable Balloon diameter¹

Motiva Implants Ergonomix2[®] Round Matrix


MINI



DEMI



FULL



Motiva[®] Inflatable Balloon 80
REF MIB-0080



Motiva[®] Inflatable Balloon 90
REF MIB-0090



Motiva[®] Inflatable Balloon 100
REF MIB-00100



INFLATABLE BALLOON	⊙ (cm) ←A→	Catalogue	D (cm) ←B→	V (cc)	Catalogue	D (cm) ←B→	V (cc)	Catalogue	D (cm) ←B→	V (cc)
80 mm	8.00	E2SM-95	2.1	95	E2SD-125	3.0	125	E2SF-135	3.3	135
	8.50	E2SM-110	2.2	110	E2SD-150	3.1	150	E2SF-155	3.5	155
	9.00	E2SM-130	2.3	130	E2SD-170	3.3	170	E2SF-190	3.7	190
90 mm	9.50	E2SM-155	2.4	155	E2SD-195	3.4	195	E2SF-215	3.9	215
	9.75	E2SM-165	2.4	165	E2SD-210	3.4	210	E2SF-240	4.0	240
	10.00	E2SM-175	2.5	175	E2SD-225	3.5	225	E2SF-250	4.1	250
100 mm	10.25	E2SM-180	2.5	180	E2SD-235	3.5	235	E2SF-275	4.2	275
	10.50	E2SM-200	2.6	200	E2SD-250	3.6	250	E2SF-290	4.3	290
	10.75	E2SM-220	2.6	220	E2SD-270	3.7	270	E2SF-320	4.4	320
	11.00	E2SM-235	2.7	235	E2SD-285	3.8	285	E2SF-330	4.5	330

¹ Motiva[®] Inflatable Balloon sizes listed are some of the possible options for the Preservé procedure. Physicians should use their best judgment and medical expertise to select the Motiva[®] Inflatable Balloon that works best for them.